

# HEPES Medium Recipe

## Directions

Modification of the original recipe that changed the buffer of volvox medium to accommodate xenic cultures of LB 1662 *Glaucosphaera vacuolata*. Alkaline medium suitable for axenic and xenic cultures, especially LB 2411 *Euglena* sp. and LB 1333 *Pithophora* sp.

For 1 L Total

1. To approximately 950 mL of dH<sub>2</sub>O, add each of the components in the order specified (except vitamins) while stirring continuously.
  2. Adjust the pH to 8.2.
  3. Bring the total volume to 1 L with dH<sub>2</sub>O.
- \*For 1.5% agar medium add 15 g of agar into the flask; do not mix.
4. Cover and autoclave medium.
  5. When cooled add vitamins.
- \*For agar medium add vitamins, mix, and dispense before agar solidifies.
6. Store at refrigerator temperature.

#	Component	Amount	Stock Solution Concentration	Final Concentration
1	Ca(NO <sub>3</sub> ) <sub>2</sub> ·4H <sub>2</sub> O (Sigma C 5676)	1 mL/L	11.8 g/100 mL dH <sub>2</sub> O	0.5 mM
2	MgSO <sub>4</sub> ·7H <sub>2</sub> O (Sigma 230391)	1 mL/L	4 g/100 mL dH <sub>2</sub> O	0.16 mM
3	Na <sub>2</sub> glycerophosphate·5H <sub>2</sub> O (Sigma G 6501 )	1 mL/L	5 g/100 mL dH <sub>2</sub> O	0.16 mM
4	KCl (Fisher P 217)	1 mL/L	5 g/100 mL dH <sub>2</sub> O	0.67 mM
5	HEPES buffer (Sigma H-3375)	0.94 g/L		3.9 mM
6	<a href="#">P-IV Metal Solution</a>	6 mL/L		
7	<a href="#">Vitamin B<sub>12</sub></a>	1 mL/L		
8	<a href="#">Biotin Vitamin Solution</a>	1 mL/L		

