

Bristol-NaCl Recipe

Directions

Bristol's medium modified for 812 *Chlamydomonas moewusii*.

For 1 L Total

1. Prepare Bristol Medium through step one.
2. Add NaCl then bring total volume to 1 L with dH₂O.
*For 1.5% agar medium pour 15 g of agar into the flask; do not mix.
3. Cover and autoclave medium.
4. Store at refrigerator temperature.

#	Component	Amount	Stock Solution Concentration	Final Concentration
1	Bristol Medium	1 L		
2	NaCl (Fisher S271-500)	3.6 g/L		0.062 M