Chu's Medium Recipe

Directions

Modification of Chu No. 10 (Bold and Wynne 1978). Suitable for some xenic cultures of Cladophora, Coronastrum and Nostoc.

For 1 L Total

- 1. Add 10 mL of Chu Stock Solution to 990 mL of sterile dH₂O.
- 2. Add sterile NaHCO₃ to medium.
- 3. Store at refrigerator temperature.

#	Component	Amount	Stock Solution Concentration	Final Concentration
1	Chu Stock Solution	10 mL/L		
2	NaHCO ₃ (Fisher S 233)	1 mL/L	1.26 g/90 mL	0.17 mM