

ES/4 Enriched Seawater Recipe

Directions

Suitable for some of J.A. West's tropical green isolates that require less of the ES-enrichment in order to maintain the desired morphology and rate of growth (West 1984, pers. Comm.).

1. Aseptically add 5 mL of sterile ES Enrichment Solution per liter of Pasteurized Seawater (30 ppt).
2. Store at refrigerator temperature.

#	Component	Amount	Stock Solution Concentration	Final Concentration
1	Pasteurized Seawater	1 L		
2	Enrichment Solution for Seawater Medium	5 mL/L		