

Ochromonas Medium Recipe

Directions

Suitable for LB 1297 Poterioochromonas malhamensis, LB 1298 Ochromonas danica and LB 1300 Ochromonas minuta.

For 1 L Total

1. To approximately 900 mL of dH₂O, add each of the components in the order specified while stirring continuously.
2. Bring total volume to 1 L with dH₂O.
3. Autoclave medium.
4. Store at refrigerator temperature.

#	Component	Amount	Stock Solution Concentration	Final Concentration
1	Glucose (BACTO-dextrose or Sigma D 9634)	1 g/L		5.6 mM
2	Tryptone (Sigma T 9410))	1 g/L		
3	Yeast extract (Bacto, Difco)	1 g/L		
4	Liver extract infusion	40 mL/L		