

1/2 Enriched Seawater Recipe

Directions

For marine algae that require salinities at 15 ppt, the basic Enriched Seawater Medium is combined with sterile distilled water prior to inoculating the new culture, such that 1/2 Enriched Seawater Medium combines Enriched Seawater Medium and sterile distilled water 1:1.

For 1 L Total

1. Combine components.
2. Store at refrigerator temperature.

#	Component	Amount	Stock Solution Concentration	Final Concentration
1	Enriched Seawater Medium	500 mL		
2	sterile dH ₂ O	500 mL		