

1/2 Erdschreiber Medium Recipe

Directions

For marine algae that require salinities at 15 ppt, the basic Erdschreiber Medium is combined with sterile distilled water prior to inoculating the new culture, such that 1/2 Erdschreiber Medium combines Erdschreiber Medium and sterile distilled water 1:1.

For 1 L Total

1. Combine components.
2. Store at refrigerator temperature.

#	Component	Amount	Stock Solution Concentration	Final Concentration
1	Erdschreiber's Medium	500 mL		
2	sterile dH ₂ O	500 mL		